



RETRACTOR MAKER

PRODUCT BROCHURE OF

Ricard Abdominal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Ricard Abdominal Retractor

SKU RM-RAR-SET-1

\$229 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Self-Retaining Base Frame	Chassis with curved shanks to improve visual access
2 (1 Pair)	Detachable Solid Lateral Blades	Broad, scoop-edged side panels (Option A: 50×58 mm / Option B: 60×80 mm)
1	Detachable Central Blade	Asymmetrical flat plate (Option A: 55×70 mm / Option B: 80×90 mm)
1	Special Threaded Screw Lock	Manual rotation locking gear to retain blade expansion
1	Compact Hook Handle	Integrated minimalist grip designed to save instrument space
1	Protective Storage / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Ricard Abdominal Retractor is an internationally acclaimed abdominal surgical retractor system engineered specifically to provide safe, expansive soft-tissue and visceral displacement during major open laparotomies and pelvic operations. Designed with curved chassis shanks and a three-blade geometry, this instrument holds the abdominal wall and pelvic viscera open under locked mechanical tension, eliminating assistant hand fatigue while preserving an unobstructed visual corridor into the peritoneal cavity.

Built to accommodate diverse incision depths, this Ricard Abdominal Retractor series is cataloged across two primary multi-blade size assortments: Option A (Compact Set, \$229) consists of one 55 mm × 70 mm central blade paired with a matching pair of 50 mm × 58 mm lateral blades; and Option B (Large Set, \$249) features one 80 mm × 90 mm central blade paired with a matching pair of 60 mm × 80 mm lateral blades. Handcrafted

from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the broad scoop-edged lateral panels gather and hold bulky muscular layers cleanly, while the special threaded screw lock manual rotation gear retains blade expansion rigidly. An integrated compact hook handle saves vital operative space.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Ricard Abdominal Retractor to withstand repetitive hospital autoclave cycles, providing general surgery suites, obstetrical theaters, and emergency trauma operating rooms worldwide with exceptional mechanical stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)