



RETRACTOR MAKER

PRODUCT BROCHURE OF

Golligher Self Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Golligher Self Retaining Retractor

SKU RM-GSRR-SET-1

\$299 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Self-Retaining Abdominal Frame Chassis	Multipurpose structural ring frame for deep general surgery access
2	Detachable Solid Lateral Blades	Flat-surfaced side panels (51 mm or 102 mm depth)
2	Detachable Central Blades	Slightly curved blades contoured to internal organ structures (51 mm or 102 mm)
3	Individual Sliding Screw Locks	Mechanical ratchet locking bolts fitted to each blade to reduce fatigue
2	Ergonomic Retracting Forceps Beams	Lightweight handle profiles with smooth perimeter edges
1	Protective Storage / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Golligher Self Retaining Retractor is an internationally acclaimed abdominal and pelvic surgical retractor system engineered specifically to provide expansive, multi-directional soft-tissue and visceral displacement during colorectal, gastrointestinal, and deep pelvic surgery. Designed with an articulated circular frame chassis, this instrument allows surgeons to position multiple solid blades circumferentially, holding the abdominal wall and pelvic viscera open under locked mechanical tension without requiring manual assistant retraction.

Built to clinical framework specifications, this Golligher Self Retaining Retractor complete set (\$299) features a self-retaining abdominal frame chassis designed for deep general surgery access. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the system includes detachable solid lateral blades with

specialized edges to secure bulky tissue layers, paired with detachable central blades slightly curved to contort smoothly to internal organ structures. The system can be configured with a Shallow Blade Assortment Pack (51 mm working depth) or a Deep Blade Assortment Pack (102 mm working depth), each locked by individual sliding screw locks to reduce user fatigue.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Golligher Self Retaining Retractor to withstand repetitive hospital autoclave cycles, providing colorectal surgical suites, trauma theaters, and pelvic reconstructive centers worldwide with exceptional structural stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



LinkedIn



Facebook



Instagram



X (Twitter)



YouTube



Pinterest