



RETRACTOR MAKER

PRODUCT BROCHURE OF

Esophageal Articulating Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Esophageal Articulating Retractor

SKU RM-EAR-299-1

\$299 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
5 mm Shaft	RM-EAR-5-1	\$299
10 mm Shaft	RM-EAR-10-1	\$299

PRODUCT DESCRIPTION

The Esophageal Articulating Retractor is an indispensable laparoscopic surgical retractor and dissector engineered specifically to provide safe, non-traumatic esophageal mobilization, hiatal exposure, and retrogastric passage during bariatric and anti-reflux surgery. Designed with an articulated distal tip capable of bending past 90 degrees, this instrument wraps smoothly around the gastroesophageal junction, allowing surgeons to establish a clean retroesophageal window without compressing adjacent vagal trunks or lacerating fragile esophageal musculature.

Built to clinical framework specifications, this Esophageal Articulating Retractor model (\$299) offers versatile shaft diameter options of 5 mm or 10 mm, with working lengths spanning from 330 mm to 375 mm (overall instrument length 438 mm to 487 mm). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the working end design feature is outfitted with a dedicated tip fenestration specifically engineered for passing bariatric bands, surgical tapes, or traction Penrose drains through the retroesophageal space. The tip articulation capability bends to an angle of more than 90°, governed by an ergonomic control wheel at the handle that locks orientation rigidly.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Esophageal Articulating Retractor to withstand repetitive hospital autoclave cycles, providing bariatric operating centers, foregut surgery suites, and minimally invasive clinical centers worldwide with unmatched mechanical stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)