



RETRACTOR MAKER

PRODUCT BROCHURE OF

Martin-Simplex Cheek and Lip Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Martin-Simplex Cheek and Lip Retractor

SKU RM-MSCLR-1

\$24 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
10 cm	RM-MSCLR-1-1	\$24
4 ½ inches / 11 cm	RM-MSCLR-1-2	\$29

PRODUCT DESCRIPTION

The Martin-Simplex Cheek and Lip Retractor is an essential oral, maxillofacial, and dental surgical retractor engineered specifically to provide safe, non-traumatic displacement of the cheeks and lips during intraoral procedures. Designed with an anatomical concave working end, this instrument shields delicate mucosal tissues from rotating burs and cutting tools while maintaining an expansive visual field over the alveolar ridge.

Built to maximize structural control in confined dental corridors, this Martin-Simplex Cheek and Lip Retractor is available in two standard lengths: 10 cm (\$24) for standard oral access and 4 ½ inches (11 cm, \$29) for deeper vestibule procedures. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), its anatomical architecture incorporates a narrow-necked shank integrated with a fenestrated, concave working blade that allows fluids to drain while giving surgical assistants fatigue-free manual leverage. Its flat, solid handle prevents intraoperative slippage and reduces user hand fatigue.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Martin-Simplex Cheek and Lip Retractor to withstand repetitive hospital sterilization cycles, providing dental clinics, oral surgery centers, and facial reconstructive suites worldwide with exceptional corrosion resistance, structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)