



# RETRACTOR MAKER

## PRODUCT BROCHURE OF

Paparella-Weitlaner Self-Retaining Retractor

### CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



# Paparella-Weitlaner Self-Retaining Retractor

**SKU** RM-PWSRR-5IN-1

**\$89 USD**

## AVAILABLE SIZES & PRICING

| SIZE / OPTION                                | SKU / DETAILS                                                               | PRICE (USD) |
|----------------------------------------------|-----------------------------------------------------------------------------|-------------|
| Paparella-Weitlaner Self-Retaining Retractor | 5" (12.7 cm) overall length, 2x2 prongs (13 mm), 33 mm spread, ratchet lock | \$89        |
| Protective Sleeve / Box                      | Individual transit packaging                                                | \$89        |

## PRODUCT DESCRIPTION

The Paparella-Weitlaner Self-Retaining Retractor is a specialized otologic and microsurgical retractor engineered specifically to provide stable, hands-free soft-tissue and periosteal retraction in narrow surgical corridors. Designed with an offset low-profile finger-ring chassis, this instrument holds postauricular and scalp incisions open under mechanical tension, eliminating assistant hand fatigue while preserving an open visual field under operating microscopes.

Built in a standard configuration profile optimized for unobstructed visibility, this Paparella-Weitlaner Self-Retaining Retractor features an overall tool length of 5 inches (12.7 cm, \$89). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the working ends incorporate an asymmetrical 2x2 prong array featuring a 13.0 mm prong length with a maximum spread capability of 33.0 mm total jaw displacement. The product line provides tissue engagement options including sharp prongs to deliver a firm hold to tougher tissues, or blunt prongs for atraumatic soft tissue manipulation. Its positive ratchet lock gears hold tension reliably throughout drilling.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Paparella-Weitlaner Self-Retaining Retractor to withstand repetitive hospital autoclave cycles, delivering exceptional mechanical stability, corrosion resistance, and complete clinical durability for ENT and neurosurgical suites worldwide.



# THANK YOU

We look forward to discussing your surgical retractor requirements.

## CONTACT RETRACTOR MAKER

EMAIL

**[info@retractormaker.com](mailto:info@retractormaker.com)**

WHATSAPP

**+92 327 8345101**

WEBSITES

**[retractormaker.com](http://retractormaker.com)**  
**[retractormanufacturer.com](http://retractormanufacturer.com)**

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

## FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)