



RETRACTOR MAKER

PRODUCT BROCHURE OF

US Army Retractor, Set of Two

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



US Army Retractor, Set of Two

SKU RM-USAR-SET-21-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
US Army Retractors	21 cm overall length, double-ended 38×22 mm & 41×25 mm blades	\$49
Protective Sleeve / Box	Standard transit packaging	\$49

PRODUCT DESCRIPTION

The US Army Retractor, Set of Two is an internationally acclaimed paired surgical instrument set engineered specifically to provide versatile, multi-depth wound margin and soft-tissue displacement during general, orthopedic, and plastic procedures. Designed with contrasting right-angled blades on opposite ends of a central fenestrated handle, these retractors allow surgical teams to transition effortlessly between shallow subcutaneous tissue retraction and deeper muscular elevation without exchanging instruments.

Built to provide an absolute lifetime warranty against rust and material defects, this US Army Retractor, Set of Two kit features a fixed length of 21 cm (\$49). Handcrafted from premium operating room-grade materials, the set is highly customizable across three distinct baseline material compositions selectable via the ordering configurations. The dual-ended blade architecture incorporates a Smaller Blade Footprint measuring 38 mm long × 22 mm wide paired with a Larger Blade Footprint measuring 41 mm long × 25 mm wide, both contoured with smooth atraumatic lips that prevent tissue creep and localized bruising.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the US Army Retractor, Set of Two to withstand repetitive hospital sterilization cycles, providing outpatient clinics, trauma operating suites, and surgical centers worldwide with exceptional structural rigidity, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)