



RETRACTOR MAKER

PRODUCT BROCHURE OF

Capsule Retractors

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Capsule Retractors

SKU RM-CR-ORTHO-1

\$64 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Single Prong (Sharp)	RM-CR-ORTHO-1-1	\$64
Double Prongs (Sharp or Blunt)	RM-CR-ORTHO-1-2	\$69
Triple Prongs	RM-CR-ORTHO-1-3	\$74

PRODUCT DESCRIPTION

Capsule Retractors are indispensable orthopedic surgical instruments engineered specifically to provide secure, non-slip displacement of thick joint capsules and fibrous soft-tissue envelopes during joint arthroplasty and reconstructive surgery. Designed to navigate deep articular recesses, these instruments allow surgical teams to maintain continuous counter-traction without slippage off slick cartilaginous or osseous borders.

Built for heavy joint leverage, the Capsule Retractors series features a standard 10-inch (approx. 25.4 cm) overall length that keeps assistant hands outside the operative line of sight. Surgeons can configure the working end across three blade widths (20 mm, 22 mm, and 30 mm) paired with versatile prong layouts: a sharp Single Prong (\$64), a Double Prong (\$69) in sharp or blunt profiles, and a broad Triple Prong (\$74) layout that distributes holding forces across wide capsular bands.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures Capsule Retractors to endure rigorous mechanical torque and repetitive hospital autoclaving cycles, providing orthopedic joint replacement and trauma teams with exceptional durability and clinical stability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)