



RETRACTOR MAKER

PRODUCT BROCHURE OF

Sailer Soft Tissue Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Sailer Soft Tissue Retractor

SKU RM-SSTR-15-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Small — 15 cm, 7×25 mm blade	RM-SSTR-15-1-1	\$39
Standard — 15 cm, 10×35 mm blade	RM-SSTR-15-1-2	\$44
Medium-Large — 15 cm, 12×55 mm blade	RM-SSTR-15-1-3	\$49
Large — 15 cm, 14×70 mm blade	RM-SSTR-15-1-4	\$54
X-Large — 15 cm, 16×80 mm blade	RM-SSTR-15-1-5	\$59

PRODUCT DESCRIPTION

The Sailer Soft Tissue Retractor is an essential oral, maxillofacial, and facial plastic surgical retractor engineered specifically to provide safe, non-traumatic soft-tissue displacement in tight surgical corridors. Designed with an ultra-smooth, curved working end, this instrument allows surgeons to retract mucosal margins, cheeks, and periosteal flaps cleanly away from bone cutting burrs and osteotomes without causing localized contusions or tears.

Built to withstand frequent autoclave sterilization cycles while providing a lifetime warranty against rust and defects, this Sailer Soft Tissue Retractor configuration features an overall tool length of 15 cm (\$44), optimized for maneuvering in confined spaces. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the blade geometry configuration incorporates a handheld, single-ended smooth blade configuration designed to prevent unnecessary trauma during soft tissue manipulation. Its ergonomic flat handle provides operating assistants with steady manual leverage throughout extended procedures.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Sailer Soft Tissue Retractor to withstand repetitive hospital sterilization cycles, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for oral surgery suites, dental implant clinics, and maxillofacial centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)