



RETRACTOR MAKER

PRODUCT BROCHURE OF

Yancoskie Abdominoplasty Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Yancoskie Abdominoplasty Retractor

SKU RM-YAR-75-1

\$99 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Yancoskie Abdominoplasty Retractor	7 ½" (19.0 cm) length, 152 mm wide rake, 6 sharp in-line prongs	\$99
Protective Sleeve / Box	Individual transit packaging	\$99

PRODUCT DESCRIPTION

The Yancoskie Abdominoplasty Retractor is an indispensable aesthetic and reconstructive plastic surgery instrument engineered specifically to provide broad, non-slip cutaneous flap elevation and fascial exposure during full tummy tuck, panniculectomy, and body contouring surgery. Designed with an ultra-wide rake profile, this instrument allows surgical assistants to apply uniform manual counter-traction across the entire lower abdominal apron, elevating thick dermis and adipose tissue cleanly away from the rectus sheath without localized tissue tearing.

Built to withstand high-frequency autoclave sterilization cycles without structural degradation, this Yancoskie Abdominoplasty Retractor configuration features an overall tool length of 7 ½ inches (19.0 cm, \$99). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) to ensure exceptional resistance against environmental corrosion, its active working width features an expansive 152 mm wide rake interface. The prong architecture incorporates a single-ended design integrating 6 sharp in-line prongs that distribute traction evenly across dense subcutaneous fascial planes, preventing accidental tissue release. Its robust ergonomic handle ensures non-slip manual grip control throughout multi-hour procedures.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Yancoskie Abdominoplasty Retractor to uphold strict clinical standards under high-intensity operating lights, providing cosmetic clinics, body contouring centers, and hospital reconstructive suites worldwide with unmatched mechanical stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)