



RETRACTOR MAKER

PRODUCT BROCHURE OF

Turner-Warwick Urological Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Turner-Warwick Urological Retractor

SKU RM-TWUR-549-1

\$549 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Self-Retaining Circular Base Frame	180 mm × 140 mm rigid ring chassis designed for 360-degree visibility
6	Adjustable Claw Blades with Curved Dual Tips	13 mm × 14 mm multi-directional retraction blades for fibrous tissue
6	Quick-Slide Retention Pin Locks	Slotted manual fasteners that lock tissue tension securely
2	Ergonomic Retracting Forceps Arms	Lightweight structured handle frames forged with smooth perimeters
1	Protective Storage / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Turner-Warwick Urological Retractor is an internationally acclaimed self-retaining urological surgical retractor system engineered specifically to provide expansive, multi-directional soft-tissue and wound margin displacement during open pelvic, bladder, and reconstructive urological procedures. Designed with an anatomical ring chassis, this system allows urologic surgeons to position multiple claw blades circumferentially, holding the abdominal wall, rectus muscles, and pelvic fascia open under continuous mechanical tension without requiring manual assistant holding.

Built to clinical framework specifications, this Turner-Warwick Urological Retractor complete set (\$549) features a rigid self-retaining circular base frame with an outer footprint of 180 mm × 140 mm designed to provide 360-degree urological visibility. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the system includes 6 adjustable claw blades with curved dual tips measuring 13 mm × 14 mm at the working

end, specifically engineered to securely hook fibrous tissue sheaths and fascia without slipping. The 6 quick-slide retention pin locks feature slotted manual fasteners that lock tissue tension securely into position, while ergonomic handle frames finished with smooth perimeters prevent soft tissue tearing.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Turner-Warwick Urological Retractor to withstand repetitive hospital autoclave cycles, providing urology operating theaters, reconstructive pelvic suites, and surgical centers worldwide with exceptional mechanical stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)