



RETRACTOR MAKER

PRODUCT BROCHURE OF

Wiltse-Gelpi Self-Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Wiltse-Gelpi Self-Retaining Retractor

SKU RM-WGSRR-1

\$99 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
28.0 cm Layout — Outward-curved arms, sharp points, ratchet lock	RM-WGSRR-1-1	\$99
30.5 cm Layout — Outward-curved arms, sharp points, ratchet lock	RM-WGSRR-1-2	\$109

PRODUCT DESCRIPTION

The Wiltse-Gelpi Self-Retaining Retractor is a specialized orthopedic and spinal surgical retractor engineered specifically to provide deep, hands-free muscular retraction during paraspinal and posterolateral spinal approaches (Wiltse approach). Designed with deep-reaching outward-curved arms, this instrument anchors firmly into heavy multifidus and longissimus muscle planes, maintaining a wide, static operative corridor to the facet joints and transverse processes.

Built to withstand frequent autoclave sterilization cycles while maintaining exceptional structural strength, this Wiltse-Gelpi Self-Retaining Retractor series is offered in two overall lengths: a 28 cm layout configuration (\$99) for standard adult posterior lumbar corridors, and an extended 30.5 cm layout configuration (\$109) for deep muscular or bariatric spinal sites. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the working tip profile features outward-curved arms terminating in sharp points to enable precise tissue separation and secure anchoring without intraoperative slippage. Its multi-position mechanical ratchet lock holds counter-traction reliably throughout decortication and pedicle screw placement.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Wiltse-Gelpi Self-Retaining Retractor to endure heavy mechanical loads, providing spine centers, neurosurgical suites, and orthopedic trauma operating rooms worldwide with exceptional durability, corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)