



RETRACTOR MAKER

PRODUCT BROCHURE OF

Bookwalter StrongArm System

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Bookwalter StrongArm System

SKU RM-BSA-1

\$899 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Bookwalter StrongArm Positioning Arm	Complete Articulated 3-Joint Positioning System
1	Field Post	Available in 12" (30 cm) or 18" (45 cm) lengths

PRODUCT DESCRIPTION

The Bookwalter StrongArm System is a table-mounted support arm designed to hold Bookwalter and compatible self-retaining retractor systems steady throughout a procedure, freeing the surgical team from manual retraction. The complete articulated three-joint positioning arm allows the height, reach and angle of the attached ring frame to be adjusted and then locked in position.

The system is supplied with a field post, available in either a 12" (30 cm) or an 18" (45 cm) length, so the working height can be matched to the patient position and the depth of the exposure. Pricing runs from \$899 to \$1,099 depending on the field post length chosen; confirm your preferred configuration when you request a quote.

All components are manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) and shipped in secure protective packaging. The Bookwalter StrongArm System is supplied non-sterile, is reusable and fully autoclavable, and must be cleaned and sterilized according to your facility's validated procedures before use.

Retractor Maker manufactures and supplies the Bookwalter StrongArm System directly from its facility near Sardaar Saleem Medident Hospital, Roras Road, Sialkot 51310, Pakistan.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)