



RETRACTOR MAKER

PRODUCT BROCHURE OF

Henner Endaural Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Henner Endaural Retractor

SKU RM-HER-475-1

\$189 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
4 1/2 in (11.4 cm)	RM-HER-475-1-1	\$179
4 3/4 in (12.06 cm) — Standard	RM-HER-475-1-2	\$189

PRODUCT DESCRIPTION

The Henner Endaural Retractor is an internationally recognized otologic surgical retractor engineered specifically to provide uncompromised, hands-free exposure of the external auditory canal, tympanic membrane, and middle ear structures during endaural procedures. Designed with an anatomical curvature, this retractor spreads meatal incisions cleanly without encroaching upon the microscopic line of sight or obstructing delicate micro-instruments.

Built for tight anatomical access in confined otologic corridors, this Henner Endaural Retractor configuration features an overall tool length of 4 3/4 inches (12.06 cm, \$189). Its handle layout incorporates an ergonomically contoured handle frame designed to provide a firm, positive grip grasp to reduce surgeon hand fatigue during prolonged procedures. The finely machined mechanical lock maintains constant, non-slip tissue counter-traction, ensuring unyielding stability during delicate ossiculoplasty or tympanoplasty dissections.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Henner Endaural Retractor to strict clinical standards, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for otolaryngology and lateral skull base surgical suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)