



RETRACTOR MAKER

PRODUCT BROCHURE OF

Budde Halo Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Budde Halo Retractor

SKU RM-BHR-1

\$2,499 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Halo Ring	Double-Hinged Circular Design
3	Halo Flex Arm	9 in (229 mm)
2	Budde Halo Support Bracket	Heavy-Duty Mounting Clamps
2	Budde Halo Support Rod	5/16 in Diameter
1	Budde Halo Pattie Tray	Integrated Perforated Tray
1	Budde Halo Mini Vise	Stainless Steel Clamping Block
2	Budde Curved Retractors	Tapered Brain Retractors
2	TEW Micro Retractor Blades	Stainless Steel Micro-Blades
2	Halo Adjustment Wrenches	System Locking Tools
1	Sterilization Box	Perforated Organizer Tray

PRODUCT DESCRIPTION

The Budde Halo Retractor is a specialized double-hinged cranial retractor system engineered to establish a low-profile, 360-degree stable working perimeter around cranial openings. Attaching directly to Mayfield skull clamp uprights or table side rails via heavy-duty support rods, it isolates brain retraction forces from the table to prevent accidental tissue displacement when the operative field shifts.

Engineered for precision microneurosurgery, the Budde Halo Retractor system features a double-hinged circular halo ring that can fold into a half-ring configuration for unrestricted lateral access. The ring itself acts as an ergonomic handrest to stabilize the surgeon's forearms during delicate dissections. The set comes complete with three 9-inch flexible friction arms, TEW micro blades, tapered curved brain retractors, an integrated pattie tray, a mini vise clamp, and dedicated adjustment wrenches.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin matte finish, Retractor Maker supplies the Budde Halo Retractor non-sterile and fully autoclavable in a heavy-duty perforated aluminum sterilization container, guaranteeing surgical precision and long-term durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)