



RETRACTOR MAKER

PRODUCT BROCHURE OF

PCL Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



PCL Retractor

SKU RM-PCLR-25-1

\$99 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
PCL Retractor	25 cm (10") overall length, 21 mm tip, 4 mm prong array	\$99
Protective Sleeve / Box	Individual transit packaging	\$99

PRODUCT DESCRIPTION

The PCL Retractor is a specialized orthopedic knee surgical instrument engineered specifically to provide safe, non-traumatic displacement of the posterior cruciate ligament (PCL) and posterior capsule during total knee arthroplasty. Designed to drop behind the posterior tibial margin, this retractor pushes the proximal tibia forward while simultaneously shielding the popliteal neurovascular bundle from rotating reamers and oscillating bone saws.

Built with exceptional structural strength, this PCL Retractor configuration features an overall tool length of 25 cm (10 inches, \$99). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the working tip width incorporates 21 mm wide tips paired with a 4 mm prongs array specifically engineered to seat securely against the posterior tibial cortex, preventing instrument displacement under vigorous joint distraction. Its heavy-duty ergonomic flat handle provides operating assistants with sustained, fatigue-free manual counter-traction.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the PCL Retractor as a non-sterile, re-sterilizable operating room device built to withstand repeated hospital autoclave cycles, delivering unyielding mechanical durability, high corrosion resistance, and complete clinical dependability for joint reconstruction and trauma surgical suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)