



RETRACTOR MAKER

PRODUCT BROCHURE OF

Kolbel Soft Tissue Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Kolbel Soft Tissue Retractor

SKU RM-KSTR-8IN-1

\$119 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Kolbel Soft Tissue Retractor	Dimensions / Details: 8" (20.32 cm) overall length, curved blunt-pronged blade · Quantity Included: 1	\$119
Protective Sleeve / Box	Dimensions / Details: Individual transit packaging · Quantity Included: 1	\$119

PRODUCT DESCRIPTION

The Kolbel Soft Tissue Retractor is a specialized orthopedic surgical instrument engineered specifically to provide safe, atraumatic muscle and capsular displacement during shoulder arthroplasty and reconstructive joint surgery. Designed to navigate difficult anatomical spaces smoothly, this instrument allows surgical teams to retract the deltoid muscle, conjoined tendon, and anterior capsule without tearing tissue margins or impinging on the axillary nerve.

Built for precision handling, this dedicated Kolbel Soft Tissue Retractor configuration features an overall tool length of 8 inches (20.32 cm, \$119). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), its prong structural geometry incorporates a curved blade layout armed with blunt prongs to minimize local tissue trauma. The curved shank drops the assistant's hands safely away from the joint line, while the ergonomic handle provides fatigue-free manual leverage throughout extended procedures.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Kolbel Soft Tissue Retractor to withstand repetitive hospital autoclave cycles, providing shoulder reconstruction, sports medicine, and orthopedic trauma centers worldwide with exceptional corrosion resistance, structural rigidity, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)