



RETRACTOR MAKER

PRODUCT BROCHURE OF

T Handle Retractor Angled Blade

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



T Handle Retractor Angled Blade

SKU RM-THRAB-1

\$99 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
8.0" (20.32 cm) — 110×32 mm or 120×32 mm toothed	RM-THRAB-1-1	\$99
11.0" (27.94 cm) — 110×32 mm or 120×32 mm toothed	RM-THRAB-1-2	\$109
14.0" (35.56 cm) — 110×32 mm or 120×32 mm toothed	RM-THRAB-1-3	\$139

PRODUCT DESCRIPTION

The T Handle Retractor Angled Blade is a heavy-duty orthopedic surgical retractor engineered specifically to provide powerful, non-slip soft-tissue and muscular displacement during hip reconstruction, spine fusion, and deep pelvic surgery. Designed with an ergonomic transverse T-handle, this instrument allows surgical assistants to apply heavy manual counter-traction with a two-handed or palm-braced grip, preventing hand fatigue and slippage during prolonged joint exposures.

Built to provide excellent structural durability and backed by a lifetime warranty against rust, this T Handle Retractor Angled Blade series is available across three overall lengths: 8 inches (20.32 cm, \$99), 11 inches (27.94 cm, \$109), and an extended 14 inches (35.56 cm, \$139) for deep pelvic and bariatric corridors. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), surgeons can select between two toothed blade dimensions: 110 mm × 32 mm or 120 mm × 32 mm. The toothed blade edge bites securely into dense fibrous fascia and anchors over cortical bone margins, preventing accidental instrument displacement.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the T Handle Retractor Angled Blade to withstand heavy mechanical loads and repeated hospital autoclave cycles, providing orthopedic trauma

centers, spine suites, and pelvic surgery theaters worldwide with exceptional durability, corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)