



RETRACTOR MAKER

PRODUCT BROCHURE OF

Darling Popliteal Self-Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Darling Popliteal Self-Retaining Retractor

SKU RM-DPSRR-1

\$109 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Option A — 54 mm × 25 mm Left / 76 mm × 38 mm Right — Deep popliteal artery exposure with lateral muscle displacement	RM-DPSRR-1-1	\$109
Option B — 76 mm × 38 mm Left / 54 mm × 25 mm Right — Deep popliteal artery exposure with medial muscle displacement	RM-DPSRR-1-2	\$109

PRODUCT DESCRIPTION

The Darling Popliteal Self-Retaining Retractor is an indispensable vascular and orthopedic surgical retractor engineered specifically to provide stable, hands-free exposure of the popliteal fossa and vascular trifurcation. Designed with asymmetrical blades that match the deep anatomical contours of the posterior knee space, this instrument displaces the gastrocnemius and hamstring muscles without compressing the popliteal vessels.

Built for complete procedural adaptability, this Darling Popliteal Self-Retaining Retractor features an adjustable self-retaining frame with customizable dual blade sets (\$109). Surgeons can configure the frame with either 54 mm × 25 mm Left / 76 mm × 38 mm Right or 76 mm × 38 mm Left / 54 mm × 25 mm Right options, with standard production runs also accommodating 76 mm × 25 mm matching ratios. The integrated multi-position ratchet lock maintains targeted tissue counter-tension automatically, eliminating assistant hand fatigue during delicate vascular bypass grafting.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Darling Popliteal Self-Retaining Retractor to withstand repetitive hospital sterilization cycles, delivering high corrosion resistance, unyielding structural rigidity, and complete

clinical durability for vascular and trauma suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)