



RETRACTOR MAKER

PRODUCT BROCHURE OF

Moberg Skin Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Moberg Skin Retractor

SKU RM-MSR-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
17 cm Variant	RM-MSR-1-1	\$44
18 cm Variant	RM-MSR-1-2	\$49

PRODUCT DESCRIPTION

The Moberg Skin Retractor is a specialized handheld surgical retractor engineered specifically to provide safe, atraumatic cutaneous flap elevation and soft-tissue displacement during plastic, reconstructive, and hand surgery. Designed with an ultra-smooth, narrow working blade, this instrument allows surgeons to maintain clear visualization of underlying tendons and digital nerves without crushing delicate dermal margins.

Built for refined tactile control, this Moberg Skin Retractor series is handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) and offers two distinct length and blade configurations: a 17 cm Overall Length Variant (\$44) featuring a 5.0 mm wide × 17.0 mm long working blade layout for tight digital corridors, and an 18 cm Overall Length Variant (\$49) featuring an 8.0 mm wide × 12.0 mm long working blade layout for broader flap retraction. The smoothly contoured blade edges prevent localized ischemia, while the ergonomic flat handle ensures fatigue-free manual leverage for operating assistants.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Moberg Skin Retractor to withstand repetitive hospital sterilization cycles, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for plastic surgery suites and hand trauma centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)