



RETRACTOR MAKER

PRODUCT BROCHURE OF

Delee Universal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Delee Universal Retractor

SKU RM-DUR-1

\$129 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Delee Universal Retractor	24 cm overall length, 50 × 70 mm curved blade	\$129
Protective Sleeve / Box	Individual transit packaging	\$129

PRODUCT DESCRIPTION

The Delee Universal Retractor (frequently termed the Delee bladder retractor) is a foundational obstetrical and gynecological surgical instrument engineered specifically to provide broad, atraumatic displacement of the urinary bladder and lower abdominal wall during cesarean sections and open pelvic procedures. Designed with an expansive, smooth curved blade, this retractor isolates the bladder safely away from lower uterine hysterotomy incisions.

Built to maximize operating room exposure, this Delee Universal Retractor configuration features an overall length of 24 cm (\$129) paired with an expansive 50 mm wide by 70 mm deep fixed blade profile. The rounded, polished lip of the blade cradles the bladder atraumatically, preventing mucosal bruising or puncture under heavy counter-traction. Its heavy-duty ergonomic flat handle provides surgical assistants with fatigue-free leverage, ensuring that the bladder remains depressed and shielded from scalpel passes and suction tips.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Delee Universal Retractor to withstand repetitive hospital sterilization cycles, providing labor and delivery units and gynecological operating theaters worldwide with exceptional durability and clinical reliability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)