



RETRACTOR MAKER

PRODUCT BROCHURE OF

Hohmann Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Hohmann Retractor

SKU RM-HR-12IN-1

\$74 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
15 mm Wide Blade — 12" (30.48 cm)	RM-HR-12IN-1-1	\$74
17 mm Wide Blade — 12" (30.48 cm)	RM-HR-12IN-1-2	\$79

PRODUCT DESCRIPTION

The Hohmann Retractor is an indispensable orthopedic bone leverage instrument engineered specifically to provide reliable tissue displacement and bone exposure during joint arthroplasty and fracture stabilization. Designed with a curved shank and an anatomical beak tip that hooks over bone margins, this retractor allows surgical assistants to leverage large muscle masses away from the bone without slipping into neurovascular zones.

Built to closely fit varied skeletal boundaries, this Hohmann Retractor maintains a unified longitudinal layout profile of 12 inches (30.48 cm). The instrument is manufactured with two selectable blade width profiles: 15 mm wide (\$74) for targeted bone leverage and 17 mm wide (\$79) for broader muscular displacement. Its contoured tip wraps securely behind the femoral or tibial cortex, while the heavy-duty handle provides surgical assistants with sustained pulling leverage and reduced hand fatigue.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Hohmann Retractor to endure heavy orthopaedic torque and repeated hospital autoclave cycles, providing joint replacement and trauma surgery teams worldwide with unyielding clinical durability and performance.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)