



RETRACTOR MAKER

PRODUCT BROCHURE OF

Caspar Cervical Retractor With Ratchet

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Caspar Cervical Retractor With Ratchet

SKU RM-CCRWR-1

\$289 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
16.5 cm	RM-CCRWR-1-1	\$289
22.3 cm	RM-CCRWR-1-2	\$299

PRODUCT DESCRIPTION

The Caspar Cervical Retractor With Ratchet is a premier spinal surgical instrument engineered specifically to provide stable, self-retaining soft-tissue exposure during anterior cervical spine interventions. Designed with a multi-jointed articulating frame, this retractor adapts smoothly to patient neck contours while shielding the carotid sheath, trachea, and esophagus from accidental intraoperative trauma.

Built for complete procedural adaptability, this Caspar Cervical Retractor With Ratchet platform is manufactured in two selectable overall lengths: a standard 16.5 cm model (\$289) and an extended 22.3 cm model (\$299) designed for broad operative exposure or muscular patient profiles. Its precision mechanical ratchet track allows spine surgeons to spread retractor blades in micro-increments, locking firmly under tissue counter-tension to eliminate manual assistant fatigue throughout discectomy and cage placement.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Caspar Cervical Retractor With Ratchet to withstand repetitive hospital autoclave cycles, delivering exceptional mechanical rigidity, smooth ratchet control, and complete clinical dependability for cervical spine suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)