



RETRACTOR MAKER

PRODUCT BROCHURE OF

Daniel Endoscopic Forehead Nerve Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Daniel Endoscopic Forehead Nerve Retractor

SKU RM-DEFNR-1

\$54 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Daniel Endoscopic Forehead Nerve Retractor	9 ¾" (17 cm) reach, Curved Left tip profile	\$54
Protective Sleeve / Box	Individual transit packaging	\$54

PRODUCT DESCRIPTION

The Daniel Endoscopic Forehead Nerve Retractor is a specialized aesthetic and plastic surgery instrument engineered specifically to provide safe, non-traumatic nerve mobilization and tissue retraction during endoscopic forehead and brow lift procedures. Designed to work in tandem with rigid endoscopes, this instrument allows surgeons to dissect and elevate subperiosteal tissue planes while protecting delicate supraorbital and supratrochlear neurovascular bundles.

Built for dynamic access in confined fields, this Daniel Endoscopic Forehead Nerve Retractor configuration features an overall reach of 9 ¾ inches (17 cm, \$54). Its working end incorporates a specialized Curved Left shape specifically engineered to successfully navigate and hold back local soft tissues around intricate nerve bundles without causing localized stretch neuropraxia or vascular tearing. The ergonomic flat handle provides tactile sensitivity and non-slip control during delicate endoscopic tunneling.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker supplies the Daniel Endoscopic Forehead Nerve Retractor non-sterile and fully autoclavable inside protective packaging, providing plastic and reconstructive surgeons worldwide with superior durability, balance, and procedural precision.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)