



RETRACTOR MAKER

PRODUCT BROCHURE OF

AORI Patellar Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



AORI Patellar Retractor

SKU RM-APR-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
5" (12.7 cm), 10 mm working end	RM-APR-1-1	\$69
7" (17.8 cm), 10 mm working end	RM-APR-1-2	\$79

PRODUCT DESCRIPTION

The AORI Patellar Retractor is an essential orthopedic knee instrument engineered to facilitate secure patellar subluxation or eversion during primary and revision total knee arthroplasty. By anchoring firmly around the patella, this retractor maintains complete visual clearance of the trochlear groove and patellofemoral joint space without placing excessive focal stress on the patellar tendon.

Designed to maximize operative holding stability, the AORI Patellar Retractor is available in two selectable lengths: a compact 5-inch model (\$69) and a 7-inch extended-reach model (\$79). The 10 mm working end incorporates an anatomically contoured design featuring three curved prongs positioned symmetrically on both sides of the arms, providing multi-point stabilization that prevents slipping off bony borders. Its ergonomic handle ensures comfortable manual control, reducing assistant hand fatigue during lengthy joint reconstructions.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker supplies the AORI Patellar Retractor non-sterile and fully autoclavable in protective packaging, delivering high mechanical strength and corrosion resistance for orthopedic surgery suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)