



RETRACTOR MAKER

PRODUCT BROCHURE OF

Bookwalter StrongArm GYN Surgery System

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Bookwalter StrongArm GYN Surgery System

SKU RM-BSAGYNS-1

\$1,599 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Top Mount StrongArm Plus with Rail Clamp	Articulated / Table Rail Compatible
4	Self Retaining Tilt Ratchet	Standard Multi-Angle
1	Perineal Ring with Notches	9½" × 14½" (24.1 × 36.8 cm)
2	Vaginal Lateral Blades	1" × 3" (2.5 × 7.6 cm)
1	Vaginal Posterior Blade	1" × 4½" (2.5 × 11.5 cm)
1	Deaver Blade	1" × 3½" (2.5 × 8.9 cm)
1	Deaver Blade	1" × 5" (2.5 × 12.7 cm)
1	Sterilization Box	Protective Case

PRODUCT DESCRIPTION

The Bookwalter StrongArm GYN Surgery System is a self-retaining retractor system mounted on the StrongArm Plus holder. The holder is clamped to the operating table rail and supports the retracting frame, providing stable, hands-free exposure during gynecological, vaginal and pelvic surgical procedures.

The mounting system is an articulated top mount StrongArm Plus with a table-rail-compatible clamp. It carries a perineal ring with notches measuring 9½" × 14½" (24.1 × 36.8 cm); the notched profile allows blades to be seated at defined positions around the perineal working area.

Four standard multi-angle self-retaining tilt ratchets lock to the ring and hold each blade at the chosen depth and angle. The blade complement comprises two vaginal lateral blades at 1" × 3" (2.5 × 7.6 cm), one vaginal posterior blade at 1" × 4½" (2.5 × 11.5 cm) and two Deaver blades at 1" × 3½" (2.5 × 8.9 cm) and 1" × 5" (2.5 × 12.7 cm).

All components are manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic). The Bookwalter StrongArm GYN Surgery System is brand new, reusable and fully autoclavable, and is supplied non-sterile — it must be cleaned and sterilized according to the healthcare facility's validated procedures before use. The Bookwalter StrongArm GYN Surgery System is supplied in a sterilization box that protects and organizes the components. Retractor Maker manufactures the system in Sialkot, Pakistan and supplies it directly to hospitals, surgeons and distributors, with a minimum order quantity of one full set.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)