



RETRACTOR MAKER

PRODUCT BROCHURE OF

Spring Wire Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Spring Wire Retractor

SKU RM-SWR-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
8.6 cm - 8.9 cm — 1×1, 3×3, or 4×4 prongs	RM-SWR-1-1	\$44
11.5 cm - 13.0 cm — 1×1, 3×3, or 4×4 prongs	RM-SWR-1-2	\$48
22.8 cm Macro — 1×1, 3×3, or 4×4 prongs	RM-SWR-1-3	\$60

PRODUCT DESCRIPTION

The Spring Wire Retractor is a specialized self-retaining surgical spreader engineered specifically to provide continuous, hands-free soft-tissue and wound margin exposure during microvascular, hand, plastic, and ophthalmic procedures. Designed with an elastic spring-action wire loop, this instrument exerts continuous, calibrated lateral counter-traction, freeing assistant hands while maintaining an unencumbered visual field under high surgical magnification.

Built to adapt to diverse anatomical depths, this Spring Wire Retractor series features an extensive variation matrix across six overall lengths: 8.6 cm (\$44), 8.9 cm (\$46), 11.5 cm (\$48), 12.8 cm (\$50), 13.0 cm (\$52), and 22.8 cm (\$60). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the working ends offer selectable prong layout configurations including 1×1 Prongs, 3×3 Prongs, and 4×4 Prongs arrays. Each prong tip is finely polished to anchor securely into wound margins without causing tearing or focal crush injury.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Spring Wire Retractor to withstand repetitive hospital autoclave cycles without loss of spring elasticity, providing microsurgeons and outpatient surgical theaters worldwide with exceptional corrosion resistance, structural resilience, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)