



RETRACTOR MAKER

PRODUCT BROCHURE OF

Helms Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Helms Retractor

SKU RM-HLR-11-1

\$109 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Helms Retractor	11 cm overall length, 2 smooth blunt prongs, screw lock	\$109
Protective Sleeve / Box	Individual transit packaging	\$109

PRODUCT DESCRIPTION

The Helms Retractor is a specialized microsurgical self-retaining instrument engineered specifically to provide safe, non-traumatic soft-tissue separation in narrow surgical corridors. Designed for high precision in confined fields, this retractor allows surgeons to maintain clear operative visualization during delicate ear, scalp, and facial procedures without crowding neighboring micro-instruments.

Built in a space-saving configuration profile, this Helms Retractor model features an overall tool length of 11 cm (\$109) tailored for maximum surgical control within closed operating zones. Its working end incorporates 2 smooth, atraumatic prongs designed for a gentle yet firm tissue grasp, preventing focal pressure necrosis on fragile skin margins and cartilage. The low-profile screw-lock mechanism ensures micro-incremental tension adjustments that hold securely throughout long micro-dissections.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Helms Retractor to withstand repetitive hospital autoclave cycles, providing ENT, neurosurgical, and plastic surgery suites worldwide with exceptional mechanical stability, corrosion resistance, and complete clinical durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)