



RETRACTOR MAKER

PRODUCT BROCHURE OF

Caspar Nerve Root Retractors

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Caspar Nerve Root Retractors

SKU RM-CNRR-SERIES-1

\$84 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
4 mm blade	RM-CNRR-SERIES-1-1	\$84
6 mm blade	RM-CNRR-SERIES-1-2	\$89
8 mm blade	RM-CNRR-SERIES-1-3	\$94
10 mm blade	RM-CNRR-SERIES-1-4	\$99

PRODUCT DESCRIPTION

Caspar Nerve Root Retractors are essential spinal microsurgery instruments engineered specifically to provide safe, non-traumatic mobilization and shielding of the dural sac and spinal nerve roots during posterior lumbar decompressions. Designed with an elongated 9 ½-inch (24.13 cm) reach, these instruments allow spine surgeons to work comfortably in deep interlaminar corridors without obstructing the microscopic view.

Built for delicate neural protection, the Caspar Nerve Root Retractors series features specialized downward-curved blade tips engineered to cradle nerve roots smoothly, preventing them from slipping back into the operative window or suffering injury from high-speed burs and rongeurs. The series is offered in four standardized blade widths: 4 mm (\$84), 6 mm (\$89), 8 mm (\$94), and 10 mm (\$99). The ergonomic flat handle incorporates fine tactile serrations that grant the operating surgeon exceptional micro-manipulation control.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces Caspar Nerve Root Retractors to strict clinical tolerances, ensuring high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for neurosurgical and spine centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)