



RETRACTOR MAKER

PRODUCT BROCHURE OF

Meniscal Repair Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Meniscal Repair Retractor

SKU RM-MRR-5IN-1

\$59 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Meniscal Repair Retractor	Dimensions / Details: 5" (12.7 cm) overall length, polished curved deflector blade · Quantity Included: 1	\$59
Protective Sleeve / Box	Dimensions / Details: Individual transit packaging · Quantity Included: 1	\$59

PRODUCT DESCRIPTION

The Meniscal Repair Retractor is a specialized orthopedic knee instrument engineered specifically to provide safe, non-traumatic neurovascular protection and capsule deflection during arthroscopic inside-out and outside-in meniscal repairs. Designed to seat behind the posterior knee capsule, this retractor acts as a rigid needle deflector, shielding the popliteal artery, vein, and tibial and peroneal nerves from advancing meniscal suture needles.

Built in a dedicated size configuration to optimize tactile feedback, this Meniscal Repair Retractor features an overall instrument length of 5 inches (12.7 cm, \$59). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the blade geometry incorporates a curved shank featuring an open-ended, rounded tip profile that conforms to the posterior femoral condyle and tibial plateau. Its surface profile finish is highly polished to minimize local tissue friction and allow effortless needle sliding directly out through the posteromedial or posterolateral incision.

Forged with an anti-glare satin matte handle and a mirror-polished spoon blade, Retractor Maker produces the Meniscal Repair Retractor to withstand repetitive hospital autoclave cycles, providing sports medicine specialists and knee surgeons worldwide with unyielding structural rigidity, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)