



RETRACTOR MAKER

PRODUCT BROCHURE OF

Morse Rib Spreader

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Morse Rib Spreader

SKU RM-MRS-SERIES-1

\$239 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Option A — 22 × 27 mm / 200 mm Spread	RM-MRS-SERIES-1-A	\$239
Option B — 30 × 28 mm / 210 mm Spread	RM-MRS-SERIES-1-B	\$259
Option C — 42 × 28 mm / 210 mm Spread	RM-MRS-SERIES-1-C	\$279
Option D — 55 × 28 mm / 210 mm Spread	RM-MRS-SERIES-1-D	\$299

PRODUCT DESCRIPTION

The Morse Rib Spreader is an advanced thoracic surgical retractor engineered specifically to provide versatile, non-traumatic intercostal rib separation and chest wall exposure during open thoracotomy. Designed with an ultra-lightweight parallel slide guide rail and multi-position movable blades, this instrument allows surgeons to adjust blade orientation dynamically along the rib curvature, distributing distraction forces evenly to eliminate focal bone stress and costal fractures.

Built to provide maximum clinical flexibility while minimizing total instrument weight, this Morse Rib Spreader inventory scales across four unique movable blade sizing configurations: Option A (\$239) features 22 mm × 27 mm lateral blades with a 200 mm maximum expansion track spread; Option B (\$259) features 30 mm × 28 mm blades with a 210 mm maximum track spread; Option C (\$279) features 42 mm × 28 mm blades with a 210 mm maximum track spread; and Option D (\$299) features 55 mm × 28 mm blades with a 210 mm maximum track spread. The system mounts four movable angled blades (two per shank), secured by an integrated turning screw lock that rigidifies the frame to decrease hand fatigue.

Crafted with Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) blades and high-strength anodized aluminum arms, Retractor Maker produces the Morse Rib Spreader to withstand repetitive hospital sterilization

cycles, providing cardiothoracic suites, trauma centers, and lung surgery clinics worldwide with unmatched durability, corrosion resistance, and complete autoclave dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)