



RETRACTOR MAKER

PRODUCT BROCHURE OF

Medial Ramus & Orbital Rim Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Medial Ramus & Orbital Rim Retractor

SKU RM-MRORR-29-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Medial Ramus & Orbital Rim Retractor	Dimensions / Details: 29 cm overall length, anatomical curved craniofacial blade · Quantity Included: 1	\$49
Protective Sleeve / Box	Dimensions / Details: Individual transit packaging · Quantity Included: 1	\$49

PRODUCT DESCRIPTION

The Medial Ramus & Orbital Rim Retractor is an essential maxillofacial and craniofacial surgical instrument engineered specifically to provide safe, non-traumatic soft-tissue retraction along deep bony contours. Designed with a dual-utility anatomical blade curve, this retractor allows surgeons to manipulate periorbital tissues away from the orbital rim and retract the medial pterygoid muscle cleanly from the ascending ramus during major facial reconstructions.

Built for delicate craniofacial interventions requiring extended reaching depth, this Medial Ramus & Orbital Rim Retractor configuration features an overall instrument length of 29 cm (\$49). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), its lightweight yet solid construction profile gives operating assistants superior manual leverage while reducing fatigue. The smoothly polished working blade wraps securely around cortical bone borders, preventing soft-tissue laceration and shielding delicate neurovascular bundles from rotating burs and osteotomes.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Medial Ramus & Orbital Rim Retractor to withstand repetitive hospital sterilization cycles, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for oral, maxillofacial, and plastic surgery suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)