



RETRACTOR MAKER

PRODUCT BROCHURE OF

Freeman Rake Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Freeman Rake Retractor

SKU RM-FRR-7-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
22 mm Wide Rake — 7" (18 cm)	RM-FRR-7-1-1	\$44
38 mm Wide Rake — 7" (18 cm)	RM-FRR-7-1-2	\$49

PRODUCT DESCRIPTION

The Freeman Rake Retractor is an internationally acclaimed aesthetic and reconstructive plastic surgery retractor engineered specifically to provide firm, non-slip skin flap elevation during rhytidectomy and soft-tissue procedures. Designed with multiple sharp rake prongs, this instrument allows surgeons to apply uniform manual counter-traction across cutaneous and subcutaneous planes without causing focal dermal tearing.

Built to match different incision boundaries, this standard Freeman Rake Retractor configuration maintains a total length profile of 7 inches (18 cm). The instrument is manufactured in two selectable rake widths: a 22 mm wide model (\$44) for localized facial flaps and an expansive 38 mm wide model (\$49) for broad cheek and neck undermining. Its curved shank keeps assistant hands clear of the dissection line, while the flat ergonomic handle ensures fatigue-free manual leverage throughout multi-hour procedures.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Freeman Rake Retractor to withstand repetitive hospital sterilization cycles, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for cosmetic and plastic surgery centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



LinkedIn



Facebook



Instagram



X (Twitter)



YouTube



Pinterest