



RETRACTOR MAKER

PRODUCT BROCHURE OF

Holsher Nerve Root Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Holsher Nerve Root Retractor

SKU RM-HNRR-235-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Holsher Nerve Root Retractor	23.5 cm overall length, wide descendent blunt blade	\$44
Protective Sleeve / Box	Individual transit packaging	\$44

PRODUCT DESCRIPTION

The Holsher Nerve Root Retractor is a specialized neurosurgical and spinal retractor engineered specifically to provide safe, non-traumatic mobilization and shielding of delicate nerve roots and the thecal sac during posterior spinal decompressions. Designed with an elongated slender shaft, this instrument gives spine surgeons deep access into the interlaminar space while keeping assistant hands outside the microscopic field of view.

Built for precision handling during delicate discectomies, this dedicated Holsher Nerve Root Retractor configuration features an overall tool length of 23.5 cm (\$44). Its working end incorporates a wide blade layout featuring descendent, blunt edges specifically engineered to cradle nerve roots smoothly, preventing them from slipping back into the operative window or sustaining injury from high-speed burrs and rongeurs. The wide blade surface distributes counter-traction evenly to prevent focal ischemic trauma.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker supplies the Holsher Nerve Root Retractor non-sterile and fully autoclavable inside protective packaging, providing spine surgeons and neurosurgical teams worldwide with exceptional durability, corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)