



RETRACTOR MAKER

PRODUCT BROCHURE OF

Sims Vaginal Retractor — Hollow Grip Handle

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Sims Vaginal Retractor — Hollow Grip Handle

SKU RM-SVRHGH-23-1

\$89 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
70 mm × 25 mm — 23 cm	RM-SVRHGH-23-1-1	\$89
76 mm × 38 mm — 23 cm	RM-SVRHGH-23-1-2	\$94
90 mm × 38 mm — 23 cm	RM-SVRHGH-23-1-3	\$99

PRODUCT DESCRIPTION

The Sims Vaginal Retractor — Hollow Grip Handle represents an advanced ergonomic single-ended evolution of the classic speculum retractor, engineered specifically to provide deep vaginal wall displacement while minimizing assistant hand fatigue. Designed with an elongated right-angled duckbill blade terminating in a soft retaining lip, this instrument depresses the posterior vaginal wall smoothly, establishing a wide, static operative corridor to the cervix and apical vault.

Built to offer a smooth, hygienic surface profile and an absolute lifetime warranty against rust, this Sims Vaginal Retractor — Hollow Grip Handle features an extended overall tool length of 23 cm that keeps assistant hands well outside the operative field. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a highly polished finish, the single-ended, atraumatic blunt blade profile is available in three distinct clinical sizing metrics: 70 mm × 25 mm (\$89), 76 mm × 38 mm (\$94), and 90 mm × 38 mm (\$99). Its hermetically sealed hollow grip handle drastically reduces total instrument weight while providing exceptional manual grip stability.

Forged with an anti-glare satin matte finish on handling surfaces, Retractor Maker manufactures the Sims Vaginal Retractor — Hollow Grip Handle to withstand repetitive hospital autoclave cycles, providing obstetrical, gynecological, and urogynecological surgical centers worldwide with superior durability, corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)