



RETRACTOR MAKER

PRODUCT BROCHURE OF

Bookwalter StrongArm General Surgery Retractor System

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Bookwalter StrongArm General Surgery Retractor System

SKU RM-BSAGSRS-1

\$2,199 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Top Mount StrongArm PLUS Holder	Height-Adjustable / Articulated
1	Heavy Duty Rail Clamp	Table Rail Compatible
2	Ratchet Mechanisms	Standard Non-Tilt
6	Tilt Ratchet Mechanisms	Multi-Angle
1	Segmented Ring	10½" × 14½" (26.7 × 36.8 cm)
1	Segmented Ring	10½" × 19" (26.7 × 48.3 cm)
1	Balfour Blade	2⅝" × 3" (6.7 × 7.6 cm)
2	Kelly Blades	2" × 3" (5.1 × 7.6 cm)
1	Kelly Blade	2" × 4" (5.1 × 10.2 cm)
1	Kelly Blade	2" × 5" (5.1 × 12.7 cm)
1	Kelly Blade	2" × 6" (5.1 × 15.2 cm)
1	Harrington Blade	1" × 5⅝" × 2½" (2.5 × 14.3 × 6.4 cm)
1	Gelpi Blade	Standard
1	Malleable Blade	3" × 6" (7.6 × 15.2 cm)

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
2	Malleable Blades	2" × 6" (5.1 × 15.2 cm)
1	Malleable Blade	1½" × 6" (3.8 × 15.2 cm)
1	Sterilization Box	Protective Case

PRODUCT DESCRIPTION

The Bookwalter StrongArm General Surgery Retractor System is a self-retaining retractor system mounted on the StrongArm PLUS Holder. The holder anchors to the operating table and supports the retracting frame, providing stable, hands-free surgical exposure during open general surgical procedures.

The mounting system consists of a height-adjustable, articulated top mount StrongArm PLUS Holder and a heavy-duty rail clamp compatible with the table rail. Two segmented rings are supplied — 10½" × 14½" (26.7 × 36.8 cm) and 10½" × 19" (26.7 × 48.3 cm) — so the frame can be sized and reconfigured to suit the length of the incision.

Two standard non-tilt ratchet mechanisms and six multi-angle tilt ratchets hold each blade at the required depth and angle. The blade complement includes a Balfour blade, four Kelly blades in 3", 4", 5" and 6" lengths, a Harrington blade, a Gelpi blade and four malleable blades in 3" × 6", 2" × 6" and 1½" × 6" sizes, covering a broad range of general surgical exposures.

All components are manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic). The Bookwalter StrongArm General Surgery Retractor System is brand new, reusable and fully autoclavable, and is supplied non-sterile — cleaning and sterilization must follow the healthcare facility's validated procedures before use. The Bookwalter StrongArm General Surgery Retractor System is supplied in a sterilization box for organized storage and processing. Retractor Maker manufactures the system in Sialkot, Pakistan and supplies it directly to hospitals, surgeons and distributors, with a minimum order quantity of one full set.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)