



# RETRACTOR MAKER

## PRODUCT BROCHURE OF

Kelly Retractor

### CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



# Kelly Retractor

**SKU** RM-KLR-26-1

**\$89 USD**

## AVAILABLE SIZES & PRICING

| SIZE / OPTION | SKU / DETAILS | PRICE (USD) |
|---------------|---------------|-------------|
| 55 mm × 60 mm | RM-KLR-26-1-1 | \$89        |
| 70 mm × 40 mm | RM-KLR-26-1-2 | \$94        |
| 95 mm × 60 mm | RM-KLR-26-1-3 | \$99        |

## PRODUCT DESCRIPTION

The Kelly Retractor is a classic deep-cavity surgical instrument engineered specifically to provide broad, atraumatic soft-tissue and abdominal wall displacement during major open procedures. Designed with an elongated right-angled blade terminating in a soft retaining lip, this retractor allows surgical teams to hold back heavy muscle layers, adipose tissue, and viscera cleanly away from the operative corridor.

Built for dependable deep-cavity clearance, this Kelly Retractor maintains a fixed, uniform tool length of 26 cm that keeps assistant hands well outside the surgical field. The instrument provides customizable blade surface footprints selectable from the purchase menu parameters: 55 mm × 60 mm (\$89), 70 mm × 40 mm (\$94), and an expansive 95 mm × 60 mm (\$99) profile, with general menu widths spanning from 2 to 4 cm and lengths spanning from 5 to 8 cm. Its ergonomic handle delivers comfortable manual leverage and superior stability throughout multi-hour procedures.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Kelly Retractor to withstand heavy mechanical loads and repeated hospital autoclave cycles, providing general, gynecological, and orthopedic operating theaters worldwide with unyielding clinical durability and dependability.



# THANK YOU

We look forward to discussing your surgical retractor requirements.

## CONTACT RETRACTOR MAKER

EMAIL

**[info@retractormaker.com](mailto:info@retractormaker.com)**

WHATSAPP

**+92 327 8345101**

WEBSITES

**[retractormaker.com](http://retractormaker.com)**  
**[retractormanufacturer.com](http://retractormanufacturer.com)**

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

## FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)