



RETRACTOR MAKER

PRODUCT BROCHURE OF

Mayo-Collins Retractor-Set of Two

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Mayo-Collins Retractor-Set of Two

SKU RM-MCR-SET-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Standard — 5 7/8 in, double-ended blades	RM-MCR-SET-1-1	\$49

PRODUCT DESCRIPTION

The Mayo-Collins Retractor-Set of Two is an essential paired surgical instrument set engineered specifically to provide reliable, multi-depth wound margin and soft-tissue retraction during general, orthopedic, and plastic surgical procedures. Designed with contrasting right-angled blades on opposing ends, these retractors allow surgical teams to transition effortlessly between shallow subcutaneous tissue retraction and deeper fascial elevation.

Built in a standard size profile optimized for controlled tissue displacement, this Mayo-Collins Retractor-Set of Two features an overall tool length of 15.3 cm (\$49). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) to guarantee long-life performance and corrosion resistance in active operating fields, each retractor incorporates a solid right-angled blade on one end and a fenestrated blade on the reverse end that preserves direct visual sightlines. The central flat handle provides operating assistants with steady, non-slip manual control.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Mayo-Collins Retractor-Set of Two to withstand repetitive hospital sterilization cycles, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for operating rooms, trauma units, and outpatient surgical centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)