



RETRACTOR MAKER

PRODUCT BROCHURE OF

Obwegeser Mandibular Channel Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Obwegeser Mandibular Channel Retractor

SKU RM-OMCR-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
8 mm Channel Size	RM-OMCR-1-1	\$44
10 mm Channel Size	RM-OMCR-1-2	\$49

PRODUCT DESCRIPTION

The Obwegeser Mandibular Channel Retractor is an indispensable maxillofacial surgical retractor engineered specifically to provide safe, non-traumatic soft-tissue displacement and osseous border protection during sagittal split ramus osteotomies (SSRO). Designed with an anatomical channel groove that hooks securely around the inferior cortex of the mandible, this instrument allows surgeons to perform lower border bone cuts while completely shielding the facial artery, anterior facial vein, and marginal mandibular nerve.

Built for precision intraoral clearance, this Obwegeser Mandibular Channel Retractor is available across two standardized blade width parameters: an 8 mm size (\$44) for narrow mandibular borders and a 10 mm size (\$49) for broader osseous anatomy. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) to ensure lasting strength and corrosion resistance, its smooth channel profile seats firmly against cortical bone without slipping into cervical fascial planes. The ergonomic handle delivers fatigue-free manual leverage for operating assistants.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Obwegeser Mandibular Channel Retractor to withstand repetitive hospital autoclave cycles, providing oral, maxillofacial, and facial reconstructive surgeons worldwide with exceptional mechanical stability, corrosion resistance, and complete clinical durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)