



RETRACTOR MAKER

PRODUCT BROCHURE OF

Sculpo Endoscopic Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Sculpo Endoscopic Retractor

SKU RM-SER-180-1

\$229 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Sculpo Endoscopic Retractor	180 × 25 mm blade, locking stopcock, gold-plated handle
1	Scope Channel Locking Screw	Integrates 5 mm or 10 mm endoscopes securely
1	Protective Sleeve / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Sculpo Endoscopic Retractor is a premier specialized aesthetic surgical retractor engineered specifically to provide expansive optical pocket elevation and stable scope alignment during endoscopic plastic surgery. Designed to support subcutaneous cutaneous flaps, this instrument creates a rigid optical cavity while steadying an endoscopic camera, allowing plastic surgeons to perform subgaleal and subfascial dissections with complete visual control.

Built to withstand high-frequency autoclave sterilization cycles while providing a lifetime warranty against rust and manufacturing defects, this Sculpo Endoscopic Retractor configuration features blade dimensions of 180 mm long × 25 mm wide (\$229). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument incorporates an integrated telescope support featuring a built-in locking stopcock mechanism to align and steady a 10 mm or 5 mm optical system seamlessly. The handle assembly integrates a single-ended blade chassis into a specialized gold-plated grip framework that ensures steady, fatigue-free manual leverage.

Forged with an anti-glare satin matte finish on the blade and a non-slip gold-accented handle, Retractor Maker produces the Sculpo Endoscopic Retractor to strict clinical tolerances, delivering exceptional mechanical stability, optical clarity, and complete autoclave durability for cosmetic, plastic, and endoscopic surgical centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)