



RETRACTOR MAKER

PRODUCT BROCHURE OF

Hohmann Retractors Angled

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Hohmann Retractors Angled

SKU RM-HRA-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
6 ½" Overall Length — 16.51 cm length, 4 mm angled tip	RM-HRA-1-1	\$69
11" Overall Length — 27.94 cm length, 7" handle, 17 mm angle bend	RM-HRA-1-2	\$79

PRODUCT DESCRIPTION

Hohmann Retractors Angled are specialized orthopedic surgical instruments engineered specifically to provide reliable bone leverage and soft-tissue displacement in confined anatomical spaces where straight retractors cause visual obstruction. Designed with a distinct angled shank and contoured tip beak, these instruments drop the assistant's hands safely away from the incision, preserving an unencumbered sightline into deep joint recesses.

Built to handle variable tissue depths, the Hohmann Retractors Angled series is available in two distinct longitudinal sizes. The 6 ½-inch (16.51 cm) overall length model (\$69) features a 4 mm angled tip optimized for small-joint and extremity bone exposure. The larger 11-inch (27.94 cm) overall length model (\$79) features a 7-inch handle configured to bend at a 17 mm angle, delivering substantial leverage for adult hip, knee, and long bone procedures without slipping off cortical borders.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces Hohmann Retractors Angled to withstand heavy orthopaedic torque and repeated hospital autoclave cycles, providing trauma and joint replacement surgical teams worldwide with exceptional durability and clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)