



RETRACTOR MAKER

PRODUCT BROCHURE OF

C Retractors

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



C Retractors

SKU RM-CR-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
12 ½" (31.8 cm)	RM-CR-1-1	\$69
14" (35.6 cm)	RM-CR-1-2	\$74

PRODUCT DESCRIPTION

The C Retractors are precision-engineered orthopedic surgical instruments designed specifically to provide powerful bone leverage and deep soft-tissue displacement during complex hip and pelvic reconstructions. Shaped with a distinctive curved shank and angled blade profile, these instruments allow surgical assistants to leverage large muscle masses away from bony margins while keeping the operative corridor clear of visual obstruction.

Built for demanding joint access, the C Retractors series is offered in two selectable overall lengths: 12 ½ inches (\$69) and an extended 14 inches (\$74) for bariatric or deep pelvic corridors. The working end features a 22 mm wide blade terminating in a carefully polished 4 mm round blunt tip. This specialized geometry allows surgeons to seat the retractor securely against the bone cortex, providing exceptional counter-traction stability without puncturing adjacent vascular structures or lacerating periosteum.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures C Retractors to withstand extreme mechanical loads and repeated hospital autoclaving cycles, providing joint replacement and orthopedic trauma surgeons with unyielding durability and tactile control.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)